## Mini Estufa EKmini

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#### 1. Setting point alarm:

(1) Setting upper-limit point for alarm: (When the difference between PV and SV is over upper-limit, the “Alarm Indicator” (ALM1) will light up.)

- Press **↔** twice. You will see **RL HL** in upper row and the setting point in lower row.
- Press “UP KEY” or “DOWN KEY” to set point, and press **SET** to enter the value. Press **SET** again to be back the PV/SV display.
- The initial value is 10.

(2) Setting lower-limit point for alarm: (When the difference between SV and PV is over lower-limit, the “Alarm Indicator” (ALM1) will light up.)

- Press **↔** three times. You will see **RL LL** in upper row and the setting point in lower row.
- Press “UP KEY” or “DOWN KEY” to set point, and press **SET** to enter the value. Press **SET** again to be back the PV/SV display.
- The initial value is 10.

#### 2. Setting temperature unit:

- Press **SET** more than 3 sec.
- Press **↔** once. You will see **tPUo** in upper row and the setting value in lower row.
- Press “UP KEY” or “DOWN KEY” to set temperature unit “C” or “F”, and then press **SET** to enter value. Press **SET** again to be back PV/SV display.
- The initial value is C.

#### 3. Setting PV shift (offset) value: (If the PV value is not correct, you can use this function to adjust the PV value).

- Press **SET** less than 3 sec.
- Press **↔** six times. You will see **tPof** in upper row and the setting value in lower row.
- Press “UP KEY” or “DOWN KEY” to set shift value, and press **SET** to enter the value. Press **SET** again to be back the PV/SV display.
- The initial value is 0.

#### 4. Setting Auto-tuning function:

- Press “UP KEY” or “DOWN KEY” to set SV value to be auto-tuning.
- Press **SET** less than 3 sec.
- You will see **At** in upper row and the setting value in lower row.
- Press “UP KEY” or “DOWN KEY” to choose “ON” to start or “OFF” to close auto-tuning function. When auto-tuning function is on, you can see the “AT” indicator blanking. Once the auto-tuning function finish, the light of “AT” will extinguish.
- The initial value is OFF.

Note: (Auto-tuning function is that PID controller can depend on the ambient air temperature to find the best way to reach the setting temperature and let the setting temperature keep stable.)

#### 5. Setting lock:

- Press **↔** four times. You will see **LoL** in upper row and the setting value in lower row.
- Press “UP KEY” or “DOWN KEY” to select locking status. **LoL1** can lock all settings and **LoL2** can lock others than SV; When “OFF” is selected, the lock function will be off. After selecting, press **SET** to enter the value. Press **SET** again to be back the PV/SV display.
- If you press **↔** and **SET** simultaneously, the “Lock” status will be released.
- The initial value is **LoL2**.

#### 6. During setting value, you may press **SET** anytime to be back PV/SV display.

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