



Tu Laboratorio, Tu Manera - El Analizador Veterinario AI Configurable

AI Vet 2

Analizador Morfológico

AI Vet 2

Analizador Morfológico

AI Vet 2

Analizador Morfológico



Detección de Parámetros de Muestra Múltiple



Sangre
46+ Parámetros



Heces
32+ Parámetros



Orina
21+ Parámetros



Fluidos
19+ Parámetros
(Opcional)

Contenido

1 Introducción Innovadora [1-3]

- ✓ Más de 100 certificados de Propiedad Intelectual
 - ✓ Más de 200 aplicaciones a pacientes
 - ✓ Certificación ISO9001
 - ✓ Certificación ISO 13485
 - ✓ Certificación CE obtenida
- Innovación en Diagnóstico Veterinario

2 AlVet 2 [4-10]

Analizador Morfológico


- Prueba de sangre: redefiniendo el estándar dorado
- Examinación de alta calidad de sedimento urinario
- Análisis fecal totalmente automático
- Examinación morfológica de fluido pleural y ascítico

3 Reporte Morfológico [11-20]

- Reporte Morfológico de Sangre
- Reporte Morfológico de Heces
- Reporte Morfológico de Orina
- Reporte Morfológico de Fluidos (Opcional)

INNOVADOR

Innovación en Diagnóstico Veterinario



AIvet 2

Método Tradicional

Pruebas avanzadas:

Multi especies y multi pruebas

1 Ampliamente aplicable



Perro



Gato



*Para obtener información detallada, por favor contacte a su distribuidor local.

2 Detección de parámetros multi muestra



Sangre
46+ Parámetros



Heces
32+ Parámetros



Orina
21+ Parámetros



Fluidos
19+ Parámetros
(Opcional)

3 Uso amigable



Tinción automática,
microscopía, clasificación,
conteo y separación



Análisis de un clic para
reportes estándar y
reportes morfológicos



Información
diagnóstica completa
basada en megadatos

10 Min


Prueba de frotis sanguíneo impulsada por Inteligencia Artificial

- NST & NSG & NSH: Distinción de inflamación y estrés/excitación
- RET & NRBC: Identificador de regeneración de anemia
- ETG: Evaluación de hemólisis intravascular
- SPH: Daño en la membrana de la célula roja
- Agglutinated PLT: Revela el estado de plaquetas más preciso
- Large PLT: Indicador de regeneración de plaquetas
- AGG: Problemas del sistema inmunológico



Mayor Conteo

- 10µl sangre, permite captura instantánea de 200,000 a 500,000 células, con reportes completos
- La óptica de nanoprecisión captura rápidamente más de 500 campos




Interpretación Más Precisa

- **Tecnología de microfluidos:** permite la formación de una sola capa de células sanguíneas, mejorando la precisión de la clasificación y conteo



Estándar de Oro **46+ Parámetros**

El analizador AIvet 2 provee resultados más precisos, empoderando a los veterinarios para diagnósticos más seguros.



Glóbulos Blancas (WBC)

- White blood cell count (WBC)
- Neutrophils (NEU)
- Neutrophil stab granulocyte (NST)
- Eosinophil (EOS)
- Neutrophil segmented granulocyte (NSG)
- Neutrophil hypersegmented granulocyte (NSH)
- Lymphocyte (LYM)
- Basophil (BAS)
- Monocyte (MON)

Plaquetas (PLT)


- Platelet count (PLT)
- Plateletcrit (PCT)
- Mean platelet volume (MPV)
- Large platelet (LPLT)
- Agglutinated platelet count (APLT)
- Platelet distribution width (PDW)

Glóbulos Rojos (RBC)

- Red blood cell count (RBC)
- Hemoglobin concentration (HGB)
- Mean red blood cell volume (MCV)
- Mean corpuscular hemoglobin concentration (MCHC)
- Red cell distribution width (RDW)
- Mean corpuscular hemoglobin (MCH)
- Hemoglobin distribution width (HDW)
- Hematocrit (HCT)
- Erythrocyte ghost (ETG)
- Reticulocyte (RET)
- Nucleated red blood cell (NRBC)
- Spherocytes (SPH)
- Agglutination RBC (AGG)

HECES

De ineficiencia repetitiva a eficiencia simple



Ahorro de tiempo y costos

Las pruebas se completan en 12 minutos, dando un valioso ahorro de tiempo a los veterinarios



Reportes mejorados y profesionales

Reporte de 32 parámetros, con imágenes y señales de fácil diagnóstico para los dueños de las mascotas

Informe de sospecha de huevo parasitario simplificado para un rediagnóstico veterinario




Más preciso

Reporte de evaluación cuantitativa de microbioma para un análisis eficiente del sistema gastrointestinal


Escaneo de 1000+/3000+ campos para mejorar los rangos de detección

Awalife Staining System: Fecal Examination Atlas

Huevos Parasitarios



Protozoa Intestinal




Microorganismos Patogénicos



Células



Función Digestiva



HECES

The unique automatic device allows detection of feces in a single run, during the veterinary visit

Huevos Parasitarios

- Ascaris egg
- Spirometra egg
- Hookworms egg
- Alaria egg
- Dipylidium caninum egg

Protozoa Intestinal

- Trichomonas
- Giardia cyst
- Giardia
- Coccidia
- Giardia trophozoite

Microorganismos Patógenos

- Campylobacter
- Spirochete
- Bacillus
- Yeast
- Helicobacter

Células

- RBC
- Epithelial cell
- WBC

Función Digestiva

- Starch granules
- Plant fibers
- Fat droplets
- Muscle fibers

ORINA

Examen de Sedimento Utinario de Alta Calidad




Menos Tiempo

- Adición de muestra en 15 segundos
- Reporte automático generado en 11 minutos, intervención manual mínima requerida



Más Preciso

- IA reconoce 4 tipos de casts, 6 tipos de cristales, células epiteliales, espermatozoides y mucosa
- Escaneo de 1000+ campos para mejorar los rangos de detección




AlVet 2


Analizador Morfológico

Awalife Staining System: Urine Sediment Staining Atlas


Cilindros



Cristales




Células



Microorganismos Patógenos



Otros



ORINA

Resultados más concretos, proveen los elementos clínicos más relevantes, acelerando los diagnósticos clínicos

Cilindros

- Hyaline cast
- Cellular cast
- Waxy cast
- Granular cast

Cristales

- Struvite
- Uric acid
- Cystine
- Calcium oxalate monohydrate
- Calcium oxalate dihydrate
- Calcium phosphate

Células

- Renal tubular epithelial cells
- Squamous epithelial cells
- Transitional epithelial cells
- Sperm

Microorganismos Patógenos

- Cocci
- Rods
- Yeast

Otros

- Lipid droplet
- Mucus

Exámenes morfológicos de fluidos pleurales (opcional) y líquido ascítico.

Fluidos: 19+ parámetros



Determine la naturaleza de los fluidos de cavidades corporales



Determine el estado de infección bacteriana



Permite diferentes diagnósticos de líquidos pleurales y ascíticos

1 Atlas del examen de ascitis bajo el sistema de tinción


| | | | |
|---|--|--|--|
| Células nucleadas Total Nucleated Cell (TNCC#) Total Inflammatory Cell (INC#) Total Granulocyte (GRL#) Neutrophils(NEU#) Lymphocytes (LYM#) | Macrophages (MAPC#) Mesothelial Cell(MCs#) Phagocytic Cell(PHC#) Unclassified Nucleated Cells(UCC#) | Eritrocitos Total Red Blood Cell Count (RBC#) Packed Cell Volume (PCV%) | Microorganismos Bacilli (BAC#) Cocci (COS#) |
| | | | |




Bacilos




Coco




Linfocito




Neutrófilo




Eosinófilo



Macrófago




Célula fagocítica




Glóbulos rojos



Grupo agregado de células



Granulocito degenerativo






Célula mesotelial

No.: _____ L I S: _____ Doctor: _____ Sample: _____ Owner: _____
 Pet name: _____ Species: _____ Gender: _____ Pet age: _____ Weight: _____


Parameters

| Detection items | Result | Unit | Reference | Low | Normal | High |
|-----------------|--------|---------------------|---------------|-----|--------|------|
| 1.WBC | 125.52 | 10 ⁹ /L | 3.50-17.90 | | | |
| 1-1.NEU# | 0.79 | 10 ⁹ /L | 2.30-12.58 | | | |
| 1-2.NST# | 0.09 | 10 ⁹ /L | 0.00-0.80 | | | |
| 1-3.NSG# | 0.70 | 10 ⁹ /L | 2.30-12.50 | | | |
| 1-4.NHG(NSH#) | 0.00 | 10 ⁹ /L | 0.00-0.30 | | | |
| 1-5.LYM# | 105.71 | 10 ⁹ /L | 0.73-6.60 | | | |
| 1-6.MON# | 18.30 | 10 ⁹ /L | 0.00-0.90 | | | |
| 1-7.EOS# | 0.73 | 10 ⁹ /L | 0.00-1.20 | | | |
| 1-8.BAS# | 0.00 | 10 ⁹ /L | 0.00-0.12 | | | |
| 1-9.NEU% | 0.63 | % | 38.00-80.00 | | | |
| 1-10.NST/WBC% | 0.07 | % | 0.00-10.00 | | | |
| 1-11.NSG% | 0.55 | % | 35.00-75.00 | | | |
| 1-12.NSH/WBC% | 0.00 | % | 0.00-3.00 | | | |
| 1-13.LYM% | 84.22 | % | 20.00-50.00 | | | |
| 1-14.MON% | 14.58 | % | 1.00-7.20 | | | |
| 1-15.EOS% | 0.58 | % | 1.00-11.20 | | | |
| 1-16.BAS% | 0.00 | % | 0.00-0.20 | | | |
| 1-17.NST/NEU% | 11.54 | % | 0.00-15.00 | | | |
| 1-18.NSH/NEU% | 0.00 | % | 0.00-4.00 | | | |
| 2.RBC | 3.06 | 10 ¹² /L | 5.60-12.60 | | | |
| 2-1.HGB | 54.26 | g/L | 98.00-178.00 | | | |
| 2-2.HCT | 14.36 | % | 26.00-47.00 | | | |
| 2-3.MCV | 46.88 | fL | 38.70-52.50 | | | |
| 2-4.MCH | 17.72 | pg | 11.80-16.50 | | | |
| 2-5.MCHC | 377.92 | g/L | 280.00-380.00 | | | |
| 2-6.RDW-SD | 20.00 | fL | 16.00-31.90 | | | |
| 2-7.RDW-CV | 14.52 | % | 15.50-24.20 | | | |
| 2-8.HDW-SD | 11.00 | g/L | 5.80-9.80 | | | |
| 2-9.HDW-CV | 20.85 | % | 13.20-23.00 | | | |
| 2-10.RET# | 0.00 | 10 ⁹ /L | 0.00-9.60 | | | |
| 2-11.RET% | 0.00 | % | 0.00-0.15 | | | |
| 2-12.NRBC# | 0.00 | 10 ⁹ /L | 0.00-0.00 | | | |
| 2-13.NRBC/WBC% | 0.00 | % | 0.00-0.00 | | | |
| 2-14.ETG# | 0.04 | 10 ¹² /L | 0.00-0.06 | | | |
| 2-15.ETG% | 1.18 | % | 0.00-2.50 | | | |
| 2-16.SPH# | 0.00 | 10 ⁹ /L | 0.00-193.66 | | | |
| 2-17.SPH% | 0.00 | % | 0.00-2.71 | | | |
| 2-18.AGG# | 0.00 | 10 ⁹ /L | 0.00-0.15 | | | |
| 3.PLT | 34.42 | 10 ⁹ /L | 140.00-547.00 | | | |
| 3-1.PCT | 0.06 | % | 0.20-0.80 | | | |
| 3-2.MPV | 18.35 | fL | 8.20-21.30 | | | |
| 3-3.LPLT# | 7.46 | 10 ⁹ /L | 0.00-103.00 | | | |
| 3-4.P-LCR | 21.67 | % | 0.00-30.00 | | | |
| 3-5.APLT# | 0.20 | 10 ⁹ /L | 0.00-0.15 | | | |
| 3-6.PDW-SD | 15.00 | fL | 7.10-31.10 | | | |
| 3-7.PDW-CV | 97.90 | % | 43.70-73.40 | | | |





| | | | | |
|-----------|----------|---------|----------|---------|
| No.: | LIS: | Doctor: | Sample: | Owner: |
| Pet name: | Species: | Gender: | Pet age: | Weight: |

RBC & PLT distribution chart




HGB distribution chart






neutrophil stab granulocyte

| | | | | | |
|------|---------------|------|---------------------------|------|---------------|
| QTY: | 13 /378photos | Co.: | 0.09 x 10 ⁹ /L | PCT: | 0.07 NST#/WBC |
|------|---------------|------|---------------------------|------|---------------|



Neutrophil segmented granulocyte

| | | | | | |
|------|---------------|------|---------------------------|------|---------------|
| QTY: | 23 /378photos | Co.: | 0.70 x 10 ⁹ /L | PCT: | 0.55 NSG#/WBC |
|------|---------------|------|---------------------------|------|---------------|




Neutrophil hypersegmented granulocyte

| | | | | | |
|------|---------------|------|---------------------------|------|---------------|
| QTY: | 19 /378photos | Co.: | 0.00 x 10 ⁹ /L | PCT: | 0.00 NSH#/WBC |
|------|---------------|------|---------------------------|------|---------------|



Lymphocyte

| | | | | | |
|------|-----------------|------|-----------------------------|------|----------------|
| QTY: | 3495 /378photos | Co.: | 105.71 x 10 ⁹ /L | PCT: | 84.22 LYM#/WBC |
|------|-----------------|------|-----------------------------|------|----------------|




Monocyte

| | | | | | |
|------|----------------|------|----------------------------|------|----------------|
| QTY: | 605 /378photos | Co.: | 18.30 x 10 ⁹ /L | PCT: | 14.58 MON#/WBC |
|------|----------------|------|----------------------------|------|----------------|




Eosinophils

| | | | | | |
|------|---------------|------|---------------------------|------|---------------|
| QTY: | 24 /378photos | Co.: | 0.73 x 10 ⁹ /L | PCT: | 0.58 EOS#/WBC |
|------|---------------|------|---------------------------|------|---------------|




Basophil


| | | | | | |
|------|--------------|------|---------------------------|------|---------------|
| QTY: | 0 /378photos | Co.: | 0.00 x 10 ⁹ /L | PCT: | 0.00 BAS#/WBC |
|------|--------------|------|---------------------------|------|---------------|




Normocyte
QTY.: 4008 / 30photos | Co.: $2.30 \times 10^{12}/L$ | PCT: 78.53 Normocyte#/RBC




Normochromic
QTY.: 1997 / 48photos | Co.: $1.42 \times 10^{12}/L$ | PCT: 41.24 Normochromic#/RBC



Macrocyte
QTY.: 309 / 30photos | Co.: $0.18 \times 10^{12}/L$ | PCT: 6.05 Macrocyte#/RBC



Hyperchromic erythrocyte
QTY.: 2607 / 48photos | Co.: $1.85 \times 10^{12}/L$ | PCT: 53.84 Hyperchromic#/RBC



Microcyte
QTY.: 500 / 30photos | Co.: $0.29 \times 10^{12}/L$ | PCT: 9.80 Microcyte#/RBC



Hypochromic erythrocyte
QTY.: 18 / 48photos | Co.: $0.01 \times 10^{12}/L$ | PCT: 0.37 Hypochromic#/RBC



Reticulocyte
QTY.: 0 / 30photos | Co.: $0.00 \times 10^9/L$ | PCT: 0.00 RET#/RBC




NRBC
QTY.: 0 / 378photos | Co.: $0.00 \times 10^9/L$ | PCT: 0.00 NRBC#/RBC



Erythrocyte Ghost
QTY.: 63 / 56photos | Co.: $0.04 \times 10^{12}/L$ | PCT: 1.18 ETG#/RBC




Spherocyte
QTY.: 0 / 56photos | Co.: $0.00 \times 10^9/L$ | PCT: 0.00 SPH#/RBC



Agglutinated RBC
QTY.: 0 / 30photos | Co.: $0.00 \times 10^9/L$ | PCT: 0.00 AGG#/RBC



Large platelet
QTY.: 13 / 56photos | Co.: $7.46 \times 10^9/L$ | PCT: 21.67 LPLT#/PLT



Agglutinated platelets
QTY.: 9 / 378photos | Co.: $0.20 \times 10^9/L$ | PCT: 0.59 APLT #/(APLT #+PLT)

Diagnostic Recommendation

Single diagnosis:

- I. **[WBC>17.90]** It is common in inflammation, hematological diseases, malignant tumors and so on.
 1. **[NEU#<2.30]** It is common in serious disease consumption, poisoning, physical and chemical damage and so on.
 2. **[LYM#>6.60]** It is common in viral infection, lymphoma, lymphatic leukemia and so on.
 3. **[MON#>0.90]** It is common in chronic infectious diseases, convalescence of diseases, and the use of glucocorticoid drugs.





- II. **[RBC<5.60]** It is common in acute / chronic hemorrhagic anemia, hemolytic anemia, nutritional anemia, aplastic anemia and so on.
 1. **[HGB<98.00]** It is common in acute / chronic hemorrhagic anemia, hemolytic anemia, nutritional anemia, aplastic anemia and so on.
 2. **[HCT<26.00]** It is common in anemia or bleeding caused by various causes, and the increase in plasma volume caused by various causes.
 3. **[HDW-SD>9.80]** It is suggested that the hemoglobin content of single red blood cell is not uniform, which can be seen in hereditary red blood cell abnormality and so on.

- III. **[PLT<140.00]** It is common in sample agglutination, hemorrhage, platelet destruction, organ detention, insufficient bone marrow formation, drug induction and so on.
 1. **[PCT<0.2]** It suggests thrombocytopenia.
 2. **[APLT#>0.15]** It is common in samples where micro-agglutination is not visible to the naked eye, Commonly seen in pathological conditions such as immune-mediated thrombocytopenia, azotemia, infectious diseases, malignant tumours, heart disease, drug-induced disorders.


Combined diagnosis:

1. **[WBC>17.90, MON#>0.9]** It is suggested that chronic inflammation or the middle and later stage of inflammation.
2. **[RBC<5.60, HGB<98.00]** It is suggested that positive cell anemia is common in aplastic anemia, acute blood loss within 40 hours, hemolysis within 40 hours, low hematopoietic function, leukemia and so on.
3. **[HGB<98.00, RET#<9.60]** It is suggested that non-regenerative anemia is common in primary / secondary erythropoietic dysfunction (such as inflammation, tumor, chronic nephropathy, chronic liver disease, thyroid / adrenocortical dysfunction, etc.), iron / copper / folic acid / VB12 deficiency, lead poisoning, bone marrow fibrosis, osteosclerosis, hypoplastic anemia and so on.
4. **[LYM#>50, LYM%>70%]** It is strong suspicion of lymphoma, and it is recommended to combine imaging studies with pathological examinations for a comprehensive diagnosis.

No.: _____ LIS: _____ Doctor: _____ Sample: _____ Owner: _____
 Pet name: _____ Species: _____ Gender: _____ Pet age: _____ Weight: _____



No.: LIS: Doctor: Sample: Time:
 Pet name: Species: Gender: Pet age: Owner:



| | | | | |
|-------|-----|-------------|--------------|--|
| | | | | |
| 120um | GRA | Co.:0.24/uL | Co.:0.22/LPF | |

| | | | | | | |
|------|-----|--------------|--------------|--|--|--|
| | | | | | | |
| 40um | MAP | Co.:56.80/uL | Co.:3.34/HPF | | | |

| | | | | | |
|------|-----|-------------|--------------|--|--|
| | | | | | |
| 40um | COD | Co.:0.57/uL | Co.:0.03/HPF | | |

| | | | | | | |
|------|----|-------------|--------------|--|--|--|
| | | | | | | |
| 40um | CP | Co.:0.61/uL | Co.:0.04/HPF | | | |

| | | | | | |
|------|-----|---------------|--------------|--|--|
| | | | | | |
| 15um | COS | Co.:119.29/uL | Co.:7.03/HPF | | |

| | | | | | |
|------|-----|--------------|--------------|--|--|
| | | | | | |
| 15um | BAC | Co.:19.65/uL | Co.:1.16/HPF | | |

| | | | | | | |
|------|-----|---------------|---------------|--|--|--|
| | | | | | | |
| 15um | FAT | Co.:667.64/uL | Co.:52.42/HPF | | | |

| | | |
|-------|-----------------------|------------------|
| | | |
| 120um | Suspected cast(SCAS#) | QTY:1/960 photos |

| | | | | |
|-------------|---------------|-----------------|----------------------------------|---------|
| No.: | LIS: | Doctor: | Sample: | Owner: |
| Pet name: | Species: | Gender: | Pet age: | Weight: |
| Color: Pink | Clarity: Mild | Smell: Odorless | Protein concentration: 2.5-5g/dL | |

Parameters

| Detection items | Result | Unit | Reference | Negative | Positive |
|---|--------|---------------------|-----------|----------|----------|
| 1. Nucleated cell | | | | | |
| 1-1. Total Nucleated Cell Count (TNCC#) | 5.18 | 10 ³ /μL | 0-0 | | +++ |
| 1-2. Inflammatory Cell Count (INC#) | 0.30 | 10 ³ /μL | 0-0 | | + |
| 1-3. Total granulocyte count (GRL#) | 0.13 | 10 ³ /μL | 0-0 | | + |
| 1-4. Neutrophils (NEU#) | 0.13 | 10 ³ /μL | 0-0 | | + |
| 1-5. Degenerative neutrophil count (D-NEU#) | 0.00 | 10 ³ /μL | 0-0 | - | |
| 1-6. Neutrophils (NEU%) | 100.00 | % | | | |
| 1-7. Degenerative neutrophil count (D-NEU%) | 0.00 | % | | | |
| 1-8. Lymphocytes (LYM#) | 0.08 | 10 ³ /μL | 0-0 | | + |
| 1-9. Macrophage (MΦ#) | 0.09 | 10 ³ /μL | 0-0 | | + |
| 1-10. Granulocyte percentage (GRL#/TNCC#) | 2.49 | % | | | |
| 1-11. Lymphocytes percentage (LYM#/TNCC#) | 1.54 | % | | | |
| 1-12. Macrophage percentage (MΦ#/TNCC#) | 1.69 | % | | | |
| 1-13. Mesothelial cell count (MCs#) | 3.83 | 10 ³ /μL | 0-0 | | +++ |
| 1-14. Phagocytic cell (PHC#) | 0.00 | 10 ³ /μL | 0-0 | - | |
| 1-15. Unclassified nucleated cells (UCC#) | 4.89 | 10 ³ /μL | | | |
| 2. Erythrocytes | | | | | |
| 2-1. Red Blood Cells (RBC#) | 174.80 | 10 ³ /μL | 0-0 | | ++++ |
| 2-2. Pack Cell Volume (PCV%) | 1.14 | % | | | |
| 3. Microorganisms | | | | | |
| 3-1. Rods (BAC#) | 0.00 | /uL | 0-0 | - | |
| 3-2. Cocci (COS#) | 0.00 | /uL | 0-0 | - | |

Diagnostic Recommendation

1. [Protein 2.5-5 g/dL, TNCC#>5, TNCC#<10] Indication of modified transudate or exudate. A comprehensive evaluation should be made considering the specific protein content, specific gravity, and cellular component proportions of the sample.

- Diagnostic criteria:
- Modified transudate: Total protein between 2.5 and 5 g/dL, specific gravity between 1.017 and 1.025.
 - Exudate: Total protein > 3 g/dL, specific gravity > 1.025.

The primary difference in cellular components lies in the proportion of neutrophils, with neutrophils being the predominant cells in exudates.

Reference: "Clinical Laboratory Diagnostics for Small Animals," 5th edition.

Exudates are typically associated with inflammatory, necrotic, infectious, or malignant conditions. They can be classified as septic or non-septic, with septic exudates involving bacterial and/or fungal infections.

Septic exudates are commonly seen in cases such as purulent pleural or abdominal effusions, gastrointestinal perforations, abscesses, late-stage bile peritonitis, late-stage uroperitoneum, ruptured pyometra, pancreatitis, sepsis, and postoperative infections.


Non-septic exudates are often found in feline infectious peritonitis, tumors, physical trauma, chronic chylous effusions, eosinophilic granulomas, and similar conditions.

When bacterial or fungal infection is strongly suspected, it is recommended to perform microbial culture to confirm the diagnosis.

Modified transudates are a specialized form of transudate, characterized by an increase in protein or cellular content from non-inflammatory or mildly inflammatory sources. These are commonly associated with cardiovascular disease, tumors, feline infectious peritonitis, chylous leakage (chylothorax/abdomen), uroperitoneum, bile peritonitis, intrathoracic or intra-abdominal hemorrhage, liver dysfunction, lung lobe torsion, diaphragmatic hernia, hyperthyroidism, and glomerular nephropathy.

For accurate diagnosis, it is advised to consider the animal's life history, medical history, and physical examination findings. Key diagnostic considerations should include cellular parameters in thoracic or abdominal effusion analysis, biochemical analysis of peritoneal fluid, and urinary protein levels. Additionally, results from cardiac and abdominal imaging should be integrated for a comprehensive evaluation.

| | | | | |
|-----------|----------|---------|----------|---------|
| No.: | L I S: | Doctor: | Sample: | Owner: |
| Pet name: | Species: | Gender: | Pet age: | Weight: |



15um | Granulocytes(GRL#) | QTY: 60/288 photos | Co.: 0.13 x10³/µL

15um | Lymphocytes(LYM#) | QTY: 37/288 photos | Co.: 0.08 x10³/µL

30um | Epithelial cell(EPC) |

30um | Mesothelial cells(MEC#)(MEC#) | QTY: 1782/288 photos | Co.: 3.83 x10³/µL

15um | RBC | QTY: 21987/294 photo | Co.: 174.80 x10³/µL

15um | Unclassified nucleated cells(UCC#) | QTY: 2272/288 photos | Co.: 4.89 x10³/µL

15um | Unclassified nucleated cells(UCC#) | QTY: 2272/288 photos | Co.: 4.89 x10³/µL