

Better Vision
Better Results

COMPLETE SYSTEM

Firefly's KT2130 Veterinary Endoscopy Set is a perfect starter kit for veterinary professionals. It contains a comprehensive array of tools for a variety of clinical and surgical applications. These rugged tools were designed to work in conjunction with each other to create a complete and affordable endoscopy system.

DURABLE

Stainless steel instruments and high quality optics allow this set to handle the tough demands of a veterinary clinic.

COMPATIBLE

Firefly Wireless Endo-Camera and Wireless LED Light Source are compatible with most third-party endoscopes.

EXPANDABLE

Firefly also offers a variety of application-specific endoscopes and accessories such as flexible forceps and Wireless TV adapters. The TV adapter enables video to be displayed on any TV/Monitor wirelessly and without the use of a computer or tablet.



KT2130

Wireless Veterinary Endoscopy Set



- Complete turn-key veterinary endoscopy solution
- Ideal for ENT, Arthroscopy, Urology, Cystoscopy and many more applications
- Wirelessly transmit live video; Capture snapshots and record videos
- Flexible forceps and wireless TV adapters also available as optional accessories



Better Vision Better Results

Included in the KT2130:

ADDITIONAL INFORMATION

For more information about each of the included items, consult the following brochures:

- Firefly DE1250
- Firefly ES201
- DE1250 Accessories

WARRANTY & PRODUCT SUPPORT

All Firefly products come with a one year warranty. Product support is provided through Firefly's worldwide distribution network and also directly from its company headquarters in Massachusetts. Support is offered at no additional cost.

COMPANY BACKGROUND

Firefly has a strong history of offering high quality veterinary products to practices worldwide. First releasing their DE551 Wireless Veterinary Video Otoscope, followed by the DE1250 Wireless Endoscope Camera, available in over 40 countries, Firefly has established itself as the innovation leader in this space.

DE1250 Wireless Endo-Camera




- Completely Wireless Operation
- Captures Snapshots & Videos
- Includes Wireless Receiver, Charging Cables, and Software.

ES201 Wireless LED Light Source



- Completely Wireless Operation
- Ultra-Bright LED. 6500K White
- 4+ Hour Operation
- Includes 2 Rechargeable Batteries

ES301 Rigid Endoscope



Material: Stainless Steel
Length: 107.5mm
Diameter: 2.7mm
Angle: 0 degrees

ES301A Endoscope Jacket



Material: Stainless Steel
Length: 116mm
Diameter: 5mm at tip
Angle: 0 degrees

Specifications are subject to change in any matter and at any time without notice. The Firefly word and the Firefly logo are trademarks or trade names of Firefly Global and its affiliates in the United States and/or other countries. All rights in such names, marks or logos are reserved by Firefly Global.

Copyright © 2016 by Firefly Global, All rights reserved. Rev 1.4

www.fireflyglobal.com



Better Vision
Better Results

Firefly Veterinary Endoscopy Accessories

COMPLETE SYSTEM

Firefly's suite of stainless steel accessories for veterinary endoscopy are perfect for a variety of clinical and surgical applications. These rugged tools were designed to work in conjunction with the DE1250 Wireless Endoscope Camera, creating a complete and affordable endoscopy system.

DURABLE

Stainless steel and high quality optics allow these accessories to handle the tough demands of a veterinary clinic.

BACKGROUND

Firefly has a strong history of offering high quality veterinary products to practices worldwide. First releasing their DE551 Wireless Veterinary Video Otoscope, followed by the DE1250 Wireless Endoscope Camera, Firefly has established themselves as the innovation leader in this space, available in over 40 countries with new regions being added year after year.



ES301 Rigid Endoscope



Material = Stainless Steel
Length = 107.5mm
Diameter = 2.7mm
Angle = 0 degrees

ES301A Endoscope Jacket



Material = Stainless Steel
Length = 116mm
Diameter = 5mm at tip
Angle = 0 degrees

**Better Vision
Better Results**

ES302 Rigid Endoscope



**Material = Stainless
Length = 175mm
Diameter = 4.0mm
Angle = 0 degrees**



ES141 Flexible Alligator Forceps



**Material = Stainless Steel
Reach = 260mm
Diameter = 1.5mm**

Specifications are subject to change in any matter and at any time without notice. The Firefly word and the Firefly logo are trademarks or trade names of Firefly Global and its affiliates in the United States and/or other countries. All rights in such names, marks or logos are reserved by Firefly Global.

Copyright © 2015 by Firefly Global, All rights reserved. Rev 1.1

www.fireflyglobal.com

