

qLabs[®] Coag Panel 2

Vet PT/APTT Combo

for Professional Point of Care



Vet PT/APTT Combo

for Professional Point of Care

Feature Highlights

- Lightweight handheld meter with option of eStation docking system
- Dried-reagent strips stored at room temperature for true ready-to-use
- Excellent correlation to central laboratory results
- Versatile sample types: venous whole blood, citrate whole blood
- Smallest blood sample volume
- Stable external QC controls
- Flexible reference intervals set by operator
- Built-in barcode scanner for accuracy and efficiency
- Built-in communication port for printers or computer networks
- Patented RACE™ technology for rapid and accurate PT and APTT determination

Animal Species

- Canine (dog)
- Feline (cat)

Sample Type

- Fresh venous whole blood
- Citrate whole blood

Sample Volume

- 10 - 20 μ L

Time to Results

- < 7 minutes

Reportable Time Range

- PT: 10 to 80 seconds
- APTT: 10 to 420 seconds

Reference Intervals

- Canine (dog)
 - PT: 14 to 20 seconds
 - APTT: 94 to 123 seconds
- Feline (cat)
 - PT: 13 to 20 seconds
 - APTT: 96 to 122 seconds

Precision

- In-house lot-to-lot and day-to-day precision <5.0 CV% (mean value)

External QC Control

- Lyophilized whole blood preparations
- Level I, Level II

Operating Temperature Range

- 10° C to 35° C

Operating Humidity Range

- 10% to 90% RH

Shelf Life at Room Temperature

- 12 months at 15° C to 30° C